


STERISETS PRE-FILLED SALINE SYRINGES

Increase your efficiency






STERISETS® PRE-FILLED SALINE SYRINGES

- :: Easy to use.
- :: Clear description on the syringes.
- :: S = sterile inside and outside, to be used in a sterile field.
- :: SNS = sterile inside, non sterile outside.
- :: Barcode printed.
- :: Anti reflux.



STERISETS is well known in the medical device industry. With the experience of manufacturing pre-filled syringes for many years, Sterisets has developed a new production line for the pre-filled syringes. The whole process is performed in a cleanroom ISO 5 environment and is fully automated. All the syringes are terminally sterilized for safe usage.




COLOUR CODED CAPS AND LABELS

The colour of the closure caps and labels indicates the saline volume, helping to easily identify each size. Safe colour coded closure caps and labels :



Sterisets® saline 0,9% NaCl
pre-printed syringe label




Bar coding EAN 13 on pouch



Bar coding EAN 128 on boxes





PRE-FILLED STERILE SALINE SYRINGES

REF.	COLOUR CAP	DESCRIPTION
14361S		Syringe luer lock 5 ml, filled with 3 ml NaCl 0,9% sterile
14362S		Syringe luer Lock 5 ml NaCl 0,9% sterile
14363S		Syringe luer lock 10 ml NaCl 0,9% sterile
14364S		Syringe luer lock 20 ml NaCl 0,9% sterile

:: STERILE PACKAGING



STERIS
MEDICAL PRODUCTS

U9628 REV 0 01-09-2016



FOR FLUSHING ONLY / ALLEEN VOOR FLUSHEN / UNIQUEMENT POUR LE RINÇAGE / VAIN HUUTELUUN / ENDAST FOR SPOILING / AUSSCHLIESSLICH ZUM SPÜLEN / LEIN NA PREPLACHOVANIE / TYLKO DO PŁUKANIA / BARE SKYLING / TIKAI SKALOSANAI / ANULUT LOPUTAMISEKS / SOLO PARA IRRIGACIÓN / SOLO PER LAVAGGIO / KUN TIL GENNEMSKYLNING / PRO PROPLACHOVANÍ / SKRITA TIK PIRAPLOVIMU

Prefilled syringes 0.9% Sodium Chloride - Saline / Voorgevulde injectiespuiten 0,9% Natriumchloride
 Seringues préremplies Chlorure de Sodium 0,9% / Vorgefüllte spritze 0,9% Natriumchlorid
 Forfyllda sprutor 0,9 % Natriumklorid – Saline/ Esitäytetty ruisku 0,9 % natrium kloridi – keittosuola
 Prefylta sprøyter 0,9% Natriumklorid – Saltvann/ Predplněné injekční stříkačky 0,9% NaCl
 Strzykawka z 0,9% chlorkiem sodu – Sol fizjologiczna/ Švirkštas užpildytas 0,9% Natrio chlorido – fiziologiniu tirpalu
 Sirocs pilditas ar izotonisku 0,9% NaCl šķīdumu/ Esitildetud süstlad 0,9% Naatriumkloriid- Soolalahus
 Siringhe precaricate 0,9% sodio cloruro - Saline/ Forfylte sprøtjør 0,9% NatriumChlorid - Saline
 Jeringas precargadas con cloruro de sodio (solución salina) al 0,9% - Saline/ Predplněné injekční stříkačky 0,9% - NaCl

Sterisets Medical Products
 Zona Industrial 1 - Lote 11 a 14
 4560-164 Guilhufe, Penafiel
 Portugal www.sterisets.com

STERILE

CE 0344




:: INNER BOX

Multilingual instructions for use in each box

Easy peel closure

Colour coded label for easy recognition






PRE-FILLED STERILE NON STERILE SALINE SYRINGES

REF.	COLOUR CAP	DESCRIPTION
14361SNS		Syringe luer lock 5 ml, filled with 3 ml NaCl 0,9% sterile inside, non sterile outside
14362SNS		Syringe luer lock 5 ml NaCl 0,9% sterile inside, non sterile outside
14363SNS		Syringe luer lock 10 ml NaCl 0,9% sterile inside, non sterile outside
14364SNS		Syringe luer lock 20 ml NaCl 0,9% sterile inside, non sterile outside

:: STERILE NON STERILE PACKAGING





GB Only for flushing. Not for sterile field. FR Pour le rinçage uniquement. Ne pas utiliser dans un environnement stérile. NL Uitsluitend voor het flushen. Niet voor steriel veld. DE Nur zum Spülen. Nicht für den Sterilbereich. ES Sólo para el lavado de vías. No utilizar en campo estéril. DK Kun til skylling. Ikke til steril område. SE Endast för spolning. Inte får steril bruk. IT Solo per lavaggio. Per ambiente non sterile. CZ Pouze k proplachování. Není určeno pro sterilní pole. FI Vain huuhteluun. Ei steriili alueille. EE Ainult loputamiseks. Kasutamiseks mittesteriilses keskkonnas. LV Tikai skalošanas. Nelietot sterilā laukā. LT Tik praplovimui. Nesteriliai aplinkai. NO Kun til skylling. Ikke til sterile områder. PL Tylko do płukania. Nie stosowanie w polu jałowym. SK Iba na preplachovanie. Nevhodné do sterilného poľa.



Sterisets Medical Products
Zona Industrial 1, Lote 11 e 14
4560-164 Guilhufe, Penafiel
Portugal www.sterisets.com


CE 0344



:: INNER BOX AND CARTON



Colour coded label
for easy recognition




PRE-FILLED SALINE SYRINGES PROCESS

All the pre-filled saline syringes are easy to use, anti reflux, have clear labels and are barcode printed. They go through four steps ... Filling, Labelling, Packing and Sterilizing. The whole process is done in a cleanroom ISO 5 environment and is fully automated.



1 - FILLING



2 - LABELLING




3 - PACKING



4 - STERILIZATION

DID YOU KNOW?

Sterisets is now producing complete sets, where the pre-filled saline syringes are inserted in the pack to save even more time and money. Some examples of procedure sets: catheterization sets, dialysis on/off sets, hemodialysis sets, etc.



WOUNDCARE PRODUCTS

NON WOVENS/DRAPES/ GOWNS/SMS/PAPER

TRAYS & BOWLS

INSTRUMENTS

PROCEDURE TRAYS & SETS

MISCELLANEOUS



Sterisets International BV
P.O. Box 474
5340 AL Oss
The Netherlands

Ketelmeer 3
5347 JX Oss
The Netherlands

Phone: +31 (0) 412 667 755
E-mail: info@sterisets.eu
www.sterisets.com